

# Motorized Stage

## Controller Guidance

Controller

X

XY

Z

Horizontal plane Z

XYZ

Goniometer

Rotary

Unit


Controller

### ■Stepping motor controller



- Stepping motor controller DS102/112 ➔ P.1-197~
- Handy terminal DT100 ➔ P.1-199~
- DS102/112 Control software  
DSCONTROL-WIN ➔ P.1-200~
- Method for return to origin ➔ P.1-201~
- Cable for DS102/112 ➔ P.1-208~
- Cable for external control ➔ P.1-208~

### ■Stepping motor driver




- Stepping motor driver ➔ P.1-205~

### ■Cable




- Cable (Between stage and controller) ➔ P.1-207~
- Cable for DS102/112 ➔ P.1-208~
- Cable for external control ➔ P.1-208~
- Cable connecting diagram ➔ P.1-209~
- Motor option supplied cable ➔ P.1-211~

### ■Motor list



- Motor list ➔ P.1-213~

## Connection Example



Controller

X

XY

Z

Horizontal plane Z

XYZ


Goniometer

Rotary

Unit

Controller

## Configurations



## Stepping Motor Controller: DS102/112 Series

Instruction  
Manual

RoHS

**DS102 series**



**DS112 series**




- ① Please use our handy terminal DT100 and software DSCONTROL-WIN for programming and teaching set-up.
- ② Control the stage with electromagnetic brake is custom-made.
- ③ Not available for PG series sensor voltage 24V.
- ④ The advance verification is required when use the program for D200series.

SPEC				
Model	DS102NR DS102MS	DS102NR-IO DS102MS-IO	DS112NR DS112MS	DS112NR-IO DS112MS-IO
Dimension	140 (W) × 300 (D) × 90 (H) mm		70 (W) × 165 (D) × 155 (H) mm	
Input power	AC100~240V±10%, 50/60Hz		DC24V±10%	
Maximum power consumption	Less than 70W		Less than 2.5A	
Driver type (Divisions)	DS102NR: Normal (Full/Half) DS102MS: Micro step (1~1/250 [16 steps])		DS112NR: Normal (Full/Half) DS112MS: Micro step (1~1/250 [16 steps])	
Driver current rating		0.75A/Phase		
Number of controlled axis		2		
Coordinate setting range		±99,999,999pls		
Driving speed setting range		1~999,999pps		
Rate of rise speed setting range		1~9,999pps		
Addition-subtraction speed rate setting range		1~9,999ms		
Machine limit		2places for each axis to CW, CCW CW - Direction of CCW (Possible logic change)		
Detected proximity origin		1place for each axis(Possible logic change)		
Detected origin		1place for each axis(Possible logic change)		
How to detect origin		12 methods		
Sensor power		DC5V		
Home position		1place for each axis(Can be set up in effective area arbitrarily)		
External COM interface		RS232C : 4,800~38,400bps【D-SUB9Pin male】 USB2.0 : Full/Low Speed Only 【USB mini B connector】 Control I/O: Input 9 points(24V photo coupler), output 12 points (open collector)		
Link function		RS485(Controllable max 6-axis in DG chain)		
Programing function		8 programs (100 steps/program, Start/finish on control I/O)		
Teaching function		64 points (positioning in control I/O)		
Interpolation function		6-axis linear interpolation(between link device is easy linear interpolation)		
Universal input and output	-	Input 16 points(24V photo coupler) Output 12 points(Open collector)	-	Input 16 points(24V photo coupler) Output 12 points(Open collector)
Accessories		1 power cable and CD-ROM(USB driver)		
Weight		2.2kg		1.2kg

## DS102 panel layout

### Front panel




①Power switch

②POWER LED

③TERMINAL

\*1

### Rear panel



④RS232C

\*2

⑤USB connector

\*3

⑥LINK connector

⑦DIP switch

\*4

⑧EMS connector

\*5

⑨I/O connector for control

\*6

⑩General I/O connector

\*7

⑪Stage motor connector


\*8

⑫Fuse holder

⑬AC inlet

## DS112 panel layout

### Front panel



①Terminals for power input

②POWER LED

③TERMINAL

\*1

④RS232C connector

\*2

⑤USB connector

\*3

⑥LINK connector

⑦DIP switch

\*4

⑧EMS connector

\*5

⑨I/O connector for control

\*6

⑩General I/O connector

\*7

⑪Stage motor connector

\*8

### Specifications (Applicable for DS102/112)

\*1: DT100 Connector for DT100 handy terminal

\*2: Dsub9P male

\*3: Mini-B type

\*4: RS232C baud rate (2bit), Link (2bit),

USB ID (2bit) Set a command response (1bit)

\*5: Model type S02B—PASK-2(LF)(SN)(Manufactured by JST)

\*6: Model type 10226-52A2PL ( Manufactured by 3M)

\*7: Model type 10236-0200 ( Manufactured by 3M)

【Only DS102□□-IO, DS112□□-IO】

\*8: Model type 09-0054-00-14 (Manufactured by Binder)

\* Please refer the instruction manual about mounting from our web site.


# Motorized Stage

## Controller

- X
- XY
- Z
- Horizontal plane Z
- XYZ
- Goniometer
- Rotary
- Unit
- Controller

### Connection integration example

#### When use our recommendation stepping motor controller (DS102/112)



### Handy Terminal: DT100

RoHS

Can be operated remotely with DT100.

Can be operated continuous driving, step driving, return to origin, program execution and teaching by hand with DT100.

The display is a 16 digit × 2 line LCD.

? It may not be possible to create or edit on DT100. Available only for program number selection, setting and stop. DSCONTROL-WIN for programming is recommended.



SPEC	
<b>Model</b>	<b>DT100</b>
Display	LCD
Input power	DC24V (from main body)
Number of keys	12 keys
Cable length	1.5m
Dimensions	73 (W) × 130 (D) × 27 (H) mm
Weight	280g

## DS102/112 Control software: DSCONTROL-WIN

DSCONTROL-WIN is the DS102/112 control software used to easily set and control the DS102/DS112 stepping motor controller connected by USB or RS232C on Microsoft® Windows. Available maximum 6-axis.



\*Sample display

### Main function

- Set the parameter for each axis.
- Manual driving(continuous driving, step driving, Absolute value driving, return to origin)
- Teaching function
- Edit, upload and download of internal programs.
- Monitoring and forced output function of general I/O port

SPEC	
Model	DSCONTROL-WIN
Number of controlled axis	6-axis
Applicable interface	USB/RS232C
Applicable OS	Microsoft® Windows 2000/XP/7/8/10

\*This software cannot be operated plural start-up at the same time.

## Stage controller sample program

Required to make a program when control our stage controller from your PC.

Available free sample program from our web site.

<http://eng.surugaseiki.com/>



\*Sample display

You do not need to follow our sample program.

The final control program should be organized by yourself.

No guarantee of the motion in all environment.

# Motorized Stage

## Method for Return to Origin

Suruga's motorized stages is different from the sensor specifications depends on models. As return to origin operation is devided into types, it is necessary to choose the correct type. Selected wrong type may be operated uncorrectly. Choose your best one whatever you need according to be recommended as below.

### ■ Recommended return to origin

Recommended origin return method	Motorized Stage sensor specification			Number or stage connector pins
	Limit	ORG	NORG	
3, 4, 9, 10	—	1	—	12
5, 6, 11, 12	2	—	—	12
3, 4, 9, 10	2	1	—	12
1, 2, 7, 8	2	1	1	16, 12 (*)

\* About stages with 12-pin specifications and NORG

Type 1, 2, 7, 8 :Select the 4 sensor cable D214-2-□□A.

Type 3, 4, 9, 10 :Select the standard cable D214-2-□□.

\*Please refer to our website or catalog for electrical specifications of our stages.


### ■ The return to origin type list

Type	Motion
Type 0	No return to origin is performed
Type 1	Detect in the direction of CCW and perform detected process for CW edge(point a) of NORG signal.Next detect an edge of CCW side(point b) of ORG signal.
Type 2	Detect in the direction of CW and perform detected process for CCW edge of NORG signal.Next detect on edge of CW side (point b) of ORG signal.
Type 3	Detect in the direction of CCW and perform detected process for CCW edge(a point) of ORG signal.
Type 4	Detect in the direction of CW and perform detected process for CW edge of ORG signal.
Type 5	Detect in the direction of CCW and perform detected process for CW edge of CWLS signal.
Type 6	Detect in the direction of CW and perform detected process for CCW edge of CWLS signal.
Type 7	After finished type1, perform detected process for CCW edge of TIMING signal.
Type 8	After finished type2, perform detected process for CW edge of TIMING signal.
Type 9	After finished Type3, perform detected process for CCW edge of TIMING signal.
Type 10	After finished Type4, perform detected process for CW edge of TIMING signal.
Type 11	After finished type5, perform detected process for CCW edge of TIMING signal.
Type 12	After finished type6, perform detected process for CW edge of TIMING signal.


Return to origin sequence shows as below.

[Type 0] No return to origin is performed


[Type 1] Detect in the direction of CCW and perform detected process for CW edge(point a) of NORG signal.



**[Type 2]** Detect in the direction of CW and perform detected process for CCW edge(point a) of NORG signal.Next detect on edge of CW side (point b) of ORG signal.



**[Type 3]** Detect in the direction of CCW and perform detected process for CCW edge(a point) of ORG signal.



# Motorized Stage

## Method for Return to Origin

Controller

X

XY

Z

Horizontal plane Z

XYZ


Goniometer

Rotary


Unit

Controller


[Type 4] Detect in the direction of CW and perform detected process for CW edge of ORG signal.




[Type 5] Detect in the direction of CCW and perform detected process for CW edge(point a) of CWLS signal.




[Type 6] Detect in the direction of CW and perform detected process for CCW edge (point a) of CWLS signal.



**[Type 7]** After finished type1, perform detected process for CCW edge(point c) of TIMING signal.



**[Type 8]** After finished type2, perform detected process for CW edge(point c) of TIMING signal.



**[Type 9]** After finished Type3, perform detected process for CCW edge of TIMING signal.

**[Type 10]** After finished Type4, perform detected process for CW edge of TIMING signal.

**[Type 11]** After finished type5, perform detected process for CCW edge of TIMING signal.

**[Type 12]** After finished type6, perform detected process for CW edge of TIMING signal.

#### In case of return to origin with TIMING signal (Type7~Type12).

When the excitation condition of stepping motor is the excitation origin by stated in driver.

TIMING output becomes ON for each 7.2degree when the motor's step angle is 0.72degree.

The AND circuit that is configured with stage origin sensor and TIMING output can be detected to improving the accuracy of the origin without motor stop position variation in the origin sensor.

○	Detection start position
●	Detection finishing position
F speed	Driving speed(setting speed)
L Speed	Start speed(setting speed)
JD (JOG)	Interval detection JOG(JD=Lspeed)
LD	Stop time of limit detectio 300msec
SD	Stop time of sensor detectio 300msec

## Stepping Motor Driver

Orientalmotor

Controller

X  
XY  
Z  
Horizontal plane Z  
XYZ  
Goniometer  
Rotary  
Unit  
Controller

## ■ DC24V type input driver

Instruction Manual

RoHS



CVD507-K-A9

## What is smooth drive function?

Micro step can be driven without changing setup in controller as necessary.

\*The photo shows an image.

SPEC				
Model	CVD507-K-A9	CRD5107P	(DISCONTINUED) SD5107P3-A22 (DISCONTINUED) DFC5107P	
Compatible motor	C005C-90215P-1 PK523HPB-C15 PK523HPB-C17 PK523HPMB PK523HPMB-C1 PK525HPB PK525HPB-C1 PK525HPMB-C1 PK544PB PK544PB-C18 PK544PMB-C18 PK546PB Motor for KS451-40(Customized)	C005C-90215P-1 C7214-9015-1 C9582-9015-1 PK523HPB-C15 PK523HPB-C17 PK523HPMB PK523HPMB-C1 PK525HPB PK525HPB-C1 PK544PB-C18 PK546PB PKE545MC-A1 PKE566MC Motor for KS451-40(Customized)		
Input power	DC24V±10% 1.4A(Max)	DC24V±10% 1.4A(Max)	DC24V-36V±10% 1A(Max)	
Rated current		0.75A/Phase		
Divisions	1~1/250	Full/Half	1~1/250	
Range of operating temperature	0~50°C	0~40°C	0~40°C	
Dimension	65(W)×45(D)×21max(H)mm	65(W)×45(D)×28(H)mm	65(W)×45(D)×28(H)mm	
Weight	0.02kg	0.04kg	0.04kg	
Input signal	Photo coupler input			

## ■ Cable for DC24V input driver

SPEC			
Model	LCS04SD5	(DISCONTINUED) LCS02CFK	
Compatible driver	CVD507-K-A9/CRD5107P/SD5107P3-A22 (*End of sale)	DFC5107P (*End of sale)	
Length	600mm		

## ■ AC100V type input driver

Instruction Manual

RoHS

SPEC						
Model	(DISCONTINUED) RKD507-A					
Compatible motor	C005C-90215P-1 C7214-9015-1	C9582-9015-1 PK523HPB-C17	PK523HPMB-C1 PK525HPB	PK525HPB-C1 PK544PB-C18	PK544PMB-C18	PK546PB Motor for KS451-40(Customized)
Kinds	Standard driver					
Input power	Single phase 100V-115V ±15% 50/60Hz 1A (Max)					
Rated current	0.75A/Phase					
Divisions	1~1/250					
Range of operating temperature	0~50°C					
Dimension	45 (W) × 90 (D) × 120 (H) mm					
Weight	0.4kg					
Input signal	Photo coupler input					

## Stepping Motor Driver

**Td** 株式会社 テクノドライブ

This case type driver (Manufactured by Techno Drive Co.,Ltd.) can be connected easily.

### ■DC24V type input driver



SPEC		
Model	KR-A5MC	KR-A55MC
Compatible motor	C005C-90215P-1 PK523HPB-C15 PK523HPB-C17 PK523HPMB PK523HPMB-C1 PK525HPB PK525HPB-C1 PK525HPMB-C1 PK544PB PK544PB-C18 PK544PMB-C18 PK546PB Motor for KS451-40(Customized)	
Input power	DC20V-35V	
Rated current	0.75A/Phase	
Divisions	Full/Half	1~1/250
Range of operating temperature	0~40°C	
Dimension	93 (W) × 45 (D) × 32 (H) mm	105 (W) × 76.5 (D) × 39.5 (H) mm
Weight	0.13kg	0.22kg
Input signal	Photo coupler input	

### ■AC100V type input driver



SPEC	
Model	KR-A535M
Compatible motor	C005C-90215P-1 PK523HPB-C15 PK523HPB-C17 PK523HPMB PK523HPMB-C1 PK525HPB PK525HPB-C1 PK525HPMB-C1 PK544PB PK544PB-C18 PK544PMB-C18 PK546PB Motor for KS451-40(Customized)
Type	Standard driver
Input power	Single phase 100V-220V ±10% 50/60Hz
Rated current	0.75A/Phase
Divisions	1~1/250
Range of operating temperature	0~50°C
Dimension	42 (W) × 170 (D) × 134.3 (H) mm
Weight	0.66kg
Input signal	Photo coupler input

# Motorized Stage

## Cable

RoHS

Controller

X

XY

Z

Horizontal plane Z

XYZ

Goniometer

Rotary

Unit

Controller

The connection cable between a motorized stage and a stepping motor controller. There is the standard 2m cable and following cables.

● 2~6m cable(in 2m increments) ● One end loose cable(Loose wire on controller side) ● Robot cable

Selectable from the following cables when you purchase the motorized stage. Must be checked for cable type in code table on the page of each product.

D214-1-2E



D214-2-2E



Loose wire on controller side



Choose the connector that is either end of the connection type for using DS102/112.

### Normal cable / minimum bending radius:R33mm

Model	D214-1-2E	D214-1-4E	D214-1-6E	D214-1-2EK	D214-1-4EK	D214-1-6EK	D214-2-2E	D214-2-4E	D214-2-6E	D214-2-2EK	D214-2-4EK	D214-2-6EK		
Existence of connector	With both end connectors	One end (controller side) loose wire						With both end connectors	One end (controller side) loose wire					
Connector on controller side	09-0341-02-14 (Binder)	-						09-0341-02-14 (Binder)	-					
Stage-side connector	SRCN6A21-16S (JAE)						HR10A-10P-12S (73) (HRS)							
Cable length	2m	4m	6m	2m	4m	6m	2m	4m	6m	2m	4m	6m		
Color of insulator														

### Robot cable / Minimum bending radius:R33mm Excellent performance in bending

Model	D214-1-2R	D214-1-4R	D214-1-6R	D214-1-2RK	D214-1-4RK	D214-1-6RK	D214-2-2R	D214-2-4R	D214-2-6R	D214-2-2RK	D214-2-4RK	D214-2-6RK		
Existence of connector	With both end connectors	One end (controller side) loose wire						With both end connectors	One end (controller side) loose wire					
Connector on controller side	09-0341-02-14 (Binder)	-						09-0341-02-14 (Binder)	-					
Stage-side connector	SRCN6A21-16S (JAE)						HR10A-10P-12S (73) (HRS)							
Cable length	2m	4m	6m	2m	4m	6m	2m	4m	6m	2m	4m	6m		
Color of insulator														

### Cable for 4 sensors / Minimum bending radius:R33mm

Model	D214-2-2EA	D214-2-4EA	D214-2-6EA	D214-2-EAK	D214-2-4EAK	D214-2-6EAK	D214-2-2RA	D214-2-4RA	D214-2-6RA	D214-2-2RAK	D214-2-4RAK	D214-2-6RAK		
Existence of connector	With both end connectors	One end (controller side) loose wire						With both end connectors	One end (controller side) loose wire					
Connector on controller side	09-0341-02-14 (Binder)	-						09-0341-02-14 (Binder)	-					
Stage-side connector	HR10A-10P-12S (73) (HRS)						HR10A-10P-12S (73) (HRS)							
Cable length	2m	4m	6m	2m	4m	6m	2m	4m	6m	2m	4m	6m		
Color of insulator														

### Multi-core cables common specification


Type	Standard cable						Robot cable					
Conductor Configuration	40pieces/0.08mm [24AWG]						50pieces/0.08mm [24AWG]					
Approximate external diameter	0.65mm											
Insulator material	PVC						ETFE					
Color	Refer to the cross section											
Sheath material	PVC											
Color	Black											
Finishing external diameter	6.70±0.15mm											
Minimum bending radius	R33mm											

## DS102/112 Dedicated cable

### ●Link cable DS100-LINK□-0.5

Use for control 4 or 6-axis with linking DS102/112


Select the cable DS100-LINK2-0.5 when will be linked two of DS102/112, and DS100-LINK3-0.5 for linking three of them.



\* Link cable for 3 units.

### ●Control I/O cable DS100-CNT-2


Use this cable if DS102/112 is operated I/O control in upper controller. Loose wire on one side.



\* Refer the user's manual for support wiring color information.

### ●Universal input and output DS100-IO-2

Connect the external equipment to DS102/112 for controlling. Loose wire on one side.



\* Refer the user's manual for support wiring color information.


SPEC				
Model	DS100-LINK2-0.5	DS100-LINK3-0.5	DS100-CNT-2	DS100-IO-2
Cable length	0.5m	0.5m	2.0m	2.0m
Remarks	For linking 2 units	For linking 3 units	For control I/O (One end loose wire)	For general I/O (One end loose wire)

## External control cable

There are cables for external control: cables for USB and RS232C.

SPEC		
Model	DS100-USB-1.8	D100-R9-2
interface	USB	RS232C
Applicable controller	DS102/DS112	DS102/DS112
PC side Connector	USB USB A terminal socket	D-sub 9P Female
Cable length	1.8m	2m

### ■RS232C cable: D100-R9-2



# Motorized Stage

## Cable Connection Diagram


Controller

X  
XY  
Z  
Horizontal plane Z  
XYZ  
Goniometer  
Rotary  
Unit  
Controller

### Standard cable


#### ●D214-1-2E (K)、D214-1-4E (K)、D214-1-6E (K)

\* The one end loose cable (tailing K) has loose wire in side of 14 pin male connector (controller).



#### ●D214-2-2E (K)、D214-2-4E (K)、D214-2-6E (K)


\* The one end loose cable (tailing K) has loose wire in side of 14 pin male connector (controller).



### ■Cable for slit origin sensor(ORG2)(For PG,KX07/08/10/12,KH,KG05/07 series)

#### ●D214-2-2EA (K)、D214-2-4EA (K)、D214-2-6EA (K)


\* The one end loose cable (tailing K) has loose wire in side of 14 pin male connector (controller).



## ■ Robot cable

### ●D214-1-2R (K)、D214-1-4R (K)、D214-1-6R (K)


\* The one end loose cable (tailing K) has loose wire in side of 14 pin male connector (controller).



\*1 The terminal is open other than the electromagnetic brake.

### ●D214-2-2R (K)、D214-2-4R (K)、D214-2-6R (K)


\* The one end loose cable (tailing K) has loose wire in side of 14 pin male connector (controller).



## ■Cable for slit origin sensor(ORG2) (For PG,KX07/08/10/12,KH,KG05/07)

### ●D214-2-2RA (K)、D214-2-4RA (K)、D214-2-6RA (K)

\* The one end loose cable (tailing K) has loose wire in side of 14 pin male connector (controller).




# Motorized Stage

## Motorized Stage Optionally Supplied Cable (for Motor)

Controller

X  
XY  
Z  
Horizontal plane Z  
XYZ  
Goniometer  
Rotary  
Unit  
Controller

### ● Motor cable with electromagnetic brake (Fixed) STPO-RK2-A-3




**Discontinued**

※ The product on this page


All handling was canceled

There is no preparation of substitute items


### ● α Step motor cable



### ● Servo motor cable (Movable) SVPM-J3HF1-B-3-02S



### ● Servo motor cable (Movable) SVEM-J3HF1-B-3



## Optionally Supplied Cable αstep motor cable

### ■ CC0□OVA2F2


2 : 2m  
5 : 5m

Product Model
KS101-30L (R) PA-2A
KS101-30L (R) PA-5A

### ■ CC0□OVA2R2

2 : 2m  
3 : 3m  
5 : 5m

Product Model
KS101-30L (R) PA-2R
KS101-30L (R) PA-5R
KXC□□□□PAP
KXG06-□□□PAP
KXL06-□□□-□□PA□P
PG□□□-□□□□PAP
SXM20-T1E-D1-3
SXM26-T1E-D1-3



### ■ CC0□OVAF


2 : 2m  
5 : 5m

Product Model
KS101-30L (R) QA-2A
KS101-30L (R) QA-5A
KS102-□□□QA-2A
KS102-□□□QA-5A

### ■ CC0□OVAR


2 : 2m  
3 : 3m  
5 : 5m

Product Model
KS101-30L (R) QA-2R
KS101-30L (R) QA-5R
KS102-□□□QA-2R
KS102-□□□QA-5R
KXS18-□□□-□□QAP
SXM20-T2E-D2E-3
SXM26-T2E-D2E-3
SXM30-T2E-D2E-3
SXM30-T3E-D3E-3
SXM45-T3E-D3E-3



### ■ CC030VPF (3m)

Product Model
SXM20-T4E-D4E-3
SXM26-T4E-D4E-3
SXM30-T5E-D5E-3
SXM45-T5E-D5E-3




## Optionally Supplied Cable

With electromagnetic brake motor cable


### ■ CC030VARB (3m)

Product Model
SXM20-T2BE-D2E-3
SXM26-T2BE-D2E-3
SXM30-T2BE-D2E-3
SXM30-T3BE-D3E-3
SXM45-T3BE-D3E-3




### ■ CC030VPFB (3m)

Product Model
KXG06-□□□-MAM
KXL06□□□-□□-MA□M
PG□□□-□□□□□MAM
KXS18□□□-□□-SAM
SXM20-T4BE-D6E-3
SXM26-T4BE-D6E-3
SXM30-T4BE-D6E-3
SXM30-T5BE-D7E-3
SXM45-T5BE-D7E-3




## Motorized Stage Optionally Supplied Cable (for Sensor)


### ●CAVE-X/KXCSensor cable for motor option (Fixed) **HR10AP-S-SB-6-2**



### ●Limit sensor cable for PG series (Fixed) **PG-H-ASSY5-2000**



### ●Slit origin sensor for PG series (Micro photosensor) (Fixed) **PM-L25**



② Motor option supplied cable not available to purchase separately.

- ▶ See page P.1-037 for PG series motor
- ▶ See page P.1-051 for CAVE-X motor
- ▶ See page P.1-077 for KXL

# Motorized Stage

## Motor List

Controller

X

XY

Z

Horizontal plane Z

XYZ

Goniometer

Rotary

Unit

Controller

Motorized stages are provided with 5-lead pentagon connection stepping motor. See the table below for specifications of motor types.

### ■ No electromagnetic brake(Oriental motor)

Motor model	C005C-90215P-1	PK523HPB-C15 PK523HPB-C17	PK523HPMB PK523HPMB-C1	PK525HPB PK525HPB-C1	
Phases	5 phase				
Rated current	0.75A/Phase				
Size	□28mm	□28mm	□28mm	□28mm	
Basic step angle	0.72°	0.72°	0.36°	0.72°	
Max. Holding torque	0.048N・m	0.048N・m	0.038N・m	0.073N・m	
Roter inertial moment	$9 \times 10^{-7}$ kg・m <sup>2</sup>	$9 \times 10^{-7}$ kg・m <sup>2</sup>	$9 \times 10^{-7}$ kg・m <sup>2</sup>	$18 \times 10^{-7}$ kg・m <sup>2</sup>	
Mass	0.11kg	0.11kg	0.11kg	0.2kg	
Temperature	-10°C～+50°C				
Recommended driver	CVD507-K-A9 CRD5107P				
Stage (single axis)	KXT KXG-C KXL-C KXC-C KS101-30-C	KGB KGW KRB KRE04 KRE06 KRW04 KRW06	PK523HPB-C15(PG-C) PK523HPB-C17(KHE)	PK523HPMB(PG-E) PK523HPMB-C1 (KXG-G/KXL-G/KXC-G/ KGW-G/KRW-G)	PK525HPB(PG-D) PK525HPB-C1 (KXG-F/KXL-F/KXC06-F/KHC)

### ■ No electromagnetic brake(Oriental motor)

Motor model	PK525HPMB-C1	PK544PMB-C18	PK544PB PK544PB-C18	PK546PB
Phases	5 phase			
Rated current	0.75A/Phase			
Size	□28mm	□42mm	□42mm	□42mm
Basic step angle	0.36°	0.36°	0.72°	0.72°
Max. Holding torque	0.09N・m	0.24N・m	0.24N・m	0.42N・m
Roter inertial moment	$19 \times 10^{-7}$ kg・m <sup>2</sup>	$60 \times 10^{-7}$ kg・m <sup>2</sup>	$57 \times 10^{-7}$ kg・m <sup>2</sup>	$114 \times 10^{-7}$ kg・m <sup>2</sup>
Mass	0.2kg	0.3kg	0.3kg	0.5kg
Temperature	-10°C～+50°C			
Recommended driver	CVD507-K-A9 CRD5107P			
Stage(single axis)	KXG-H/KG-H	KX-G KS101-30-G KS102-G KG-G KS402-75G	PK544PB (KS402-180C/KRE10360) PK544PB-C18 (KS332-C/KS402-100C)	KXS-J

### ■ Motor for EtherCAT (SURUGASEIKI)

Motor model	STM28W100A
Size	□28mm
Resolution (Set to 1000P/R)	0.36°
Max. Holding torque	0.085N・m
Mass	0.12kg
Recommended driver	DS1000A-EC-28
Driver power input	DC24V±10%
Stage(single axis)	KXG-EA/KGW-EA/KRW-EA

**With electromagnetic brake (Oriental Motor Co., Ltd.)**

Motor model	PKE545MC-A1	PKE566MC	C103A-90215PM (*1)
Phase		5 Phase	
Rated current	0.35A/Phase	0.75A/Phase	0.75A/
Size	□42mm	□60mm	□28mm
Basic step angle		0.72°	
Max. Holding torque	0.27N・m	0.96N・m	0.048N・m
Rotor inertial moment	79×10 <sup>-7</sup> kg・m <sup>2</sup>	430×10 <sup>-7</sup> kg・m <sup>2</sup>	14×10 <sup>-7</sup> kg・m <sup>2</sup>
Input voltage of excitation brake	DC24V±0.5 0.08A	DC24V±0.5 0.25A	DC24V±0.5 0.08A
Excitation brake static friction torque	0.13N・m	0.48N・m	0.48N・m
Mass	0.52kg	1.2kg	0.17kg
Temperature		-10°C～+50°C	
Driver type	RKSD503M-A(C)	RKSD507M-A	—
Stage(single axis)	PG-MA KXG-MA(MB) KXL-MA	KXS-SA	KXG-MG

\*1 SURUGASEIKI original motor

**αstep (Oriental Motor Co., Ltd.)**

Motor model	ARM24SAK	AZM24AK	ARM46AC
Size	□28mm	□28mm	□42mm
Resolution (Set to 1000P/R)		0.36°	
Max. Holding torque	0.055N・m	0.095N・m	0.3N・m
Rotor inertial moment	11×10 <sup>-7</sup> kg・m <sup>2</sup>	9.2×10 <sup>-7</sup> kg・m <sup>2</sup>	58×10 <sup>-7</sup> kg・m <sup>2</sup>
Mass	0.15kg	0.15kg	0.47kg
Temperature		-10°C～+50°C	
Driver type	ARD-K	AZD-K	ARD-A
Driver power input	DC24V±10%	DC24V±5%	Single phase 100-115V 50/60Hz
Driver input current	0.9A	1.6A	2.9A
Stage(single axis)	KXC-PA KS101-30LPA(RPA) PG-PA KXG-PA KXL-PA KGW-PA/KRW-PA	KXG-ZA KGW-ZA KRW-ZA	KS101-30LQA(RQA) KS102-□LQA(RQA) KXS-QA

**AC servo motor (Mitsubishi electric co., Ltd.)**

Motor model	HG-KR053	HG-KR13
Size	□40mm	
Resolution (1Rotary)	4194304p/rev	
Speed detector	22 bits encoder	
Rated rotation speed	3000r/min	
Rated output	50W	100W
Rated torque	0.16N・m	0.32N・m
Max. Torque	0.56N・m	1.10N・m
Inertial moment	0.045×10 <sup>-4</sup> kg・m <sup>2</sup>	0.0777×10 <sup>-4</sup> kg・m <sup>2</sup>
Rated current	0.9A	0.8A
Mass	0.35kg	0.56kg
Temperature	0°C～40°C	
Driver type	MR-J4-10A	
Stage(single axis)	PG-UA/KXG-UA KXL-UA/ KGW-UA/KRW-UA	KXS-WA

**AC servo motor (Panasonic)**

Motor model	MSMF5AZL1A2
Size	□38mm
Resolution (1Rotary)	8388608p/rev
Speed detector	23bits encoder
Rated rotation speed	3000r/min
Rated output	50W
Rated torque	0.16N・m
Max. Torque	0.48N・m
Inertial moment	0.026×10 <sup>-4</sup> kg・m <sup>2</sup>
Rated current	1.1A
Mass	0.32kg
Temperature	0°C～40°C
Driver type	MADLT05SF
Stage(single axis)	KXG-UG KGW-UG KRW-UG